



Dandle PAL™

Neurodevelopmental Support Systems
Posture Alignment and a Sprinkle of Love



Small **pellet-filled bag** gives you more height, support and stability where you need it.

“**Tail**” helps to stabilize the smaller bag or any tubing or equipment by pulling it under the main part of the PAL.



“**Legs**” are extra long to provide flexion from head to toe. Longer legs give ample room to move the pellets around to individualize support for the baby.

- Use in conjunction with the Dandle ROO or Dandle WRAP when more support is needed
- Extremely versatile – use parts, pieces or multiple PALS to further individualize care
- Use two or more to help contain a baby during phototherapy or when other positioning aids are not appropriate
- Help to reinforce positioning if baby is agitated while still allowing some freedom of movement
- Provide support for tubing or equipment
- Can be used like an “extra set of hands”

Dandle PAL



SUPINE POSITION:



Two Dandle PALs are used to provide containment and boundaries entirely around the baby. This might be appropriate in conjunction with another positioning aid, like the Dandle WRAP, during phototherapy or when other positioning aids are not appropriate for medical reasons. Using the Dandle PAL to support the lower extremities will help to support a posterior pelvic tilt.



The Dandle PAL is helping to keep the head in midline, provide lateral support of the trunk, and some support for shoulder rounding. The smaller beanbag is on top of the horseshoe piece to give more height in support of the SiPAP tubing.

SIDE-LYING POSITION:



The "legs" of the Dandle PAL are stretched wide to provide a boundary for the back and head. The smaller pouch is behind the head to give a firmer surface to help prevent arching. Another PAL is used for additional support when the baby is agitated.



Again, the "legs" are stretched wide and molded behind the baby's back and head to aid in alignment of the head, neck, shoulders and trunk.

PRONE POSITION:



The Dandle PAL is used to keep the baby's legs in a tucked position. It also serves as a boundary to keep the baby from sliding down the mattress. A gel pillow is under his head and chest to keep the spine in alignment.