

Easy to use on/off
button and mode
select button



Tranquilo[®] Soothing Mat

Gentle Vibration Helps to Calm and Soothe Infants



Tranquilo Soothing Mat

- 6 vibrating modes:
 - Extra low
 - Low
 - Medium
 - High
 - Low with slow pulsing beat
 - High with fast pulsing beat
- Wipeable and reusable with soft disposable covers
- Automatic 1 hour shut off
- Size (11" x 23") fits hospital bassinets
- Battery powered for easy portability (4 AA batteries required)
- One year limited warranty
- Expected motor life of 2000 hours

Developed by a labor & delivery nurse, the Tranquilo Soothing Mat mimics mom's sounds and motions to comfort babies who are agitated, stressed or who need sensory input. The vibration levels and pulsatile rhythm choices, which imitate the intrauterine heartbeat, can easily be changed with the push of a button.

The Tranquilo Soothing Mat is an effective tool for soothing fussy and hard to console babies. In the hospital setting, clinicians find it easy to use and appreciate the different vibration settings that make the mat customizable to each baby's unique needs.

The Tranquilo Soothing Mat was featured on Shark Tank and was a semi-finalist in the 2019 Sleep Foundation Sleep Tech Awards.

Scan to see
video IFU



Ordering Information

D 20100 EA | Tranquilo Soothing Mat | Qty: 1/ea

D 20111 B25 | Tranquilo Soothing Mat Disposable Covers | Qty: 25/case

Copyright, Tranquilo Baby Holdings LLC 2020

Dandle•LION Medical | 22 Shelter Rock Lane | Danbury, CT 06810 | 203-791-9001 | www.dandleLIONmedical.com

Rev1:0820

